



Optional SmartMover XL with extended 4" platform
– available in cinnamon red only

SMARTMOVER

SPECIFICATIONS MODEL 450000

TRAVEL SURFACES

- Sidewalks, linoleum, carpeting, indoor and outdoor hard surfaces

SPEED AND POWER

- Programmable forward speed 2.3mph (3.7kph)/reverse speed 1.9mph (3.1kph)
- 3 amp charging, 24V electrical system
- (2) 33/35 amp, airline-safe, sealed batteries
- Battery capacity 40 miles (64km)*
- Onboard battery charger w/8' retractile cord
- Battery LED indicator

UNIPAC & VARI-SPEED DRIVE SYSTEM

- Maintenance-free, sealed, lifetime-lubricated, transaxle drive

MANEUVERABILITY

- 38.0-inch (96.5cm) turning radius @front wheel
- 44.0-inch (112cm) turning radius @basket corner

DIMENSIONS

- Weight w/batteries 186 lbs (84.4kgs)
- Length 49-inches (124.5cms)
- Height 37-inches (94 cms)
- Height (at top of basket) 30.5-inches (77.5 cms)
- Width 24-inches (61.0cms)
- Trimline basket dimensions:
 - Width 15.8-inches (40.1cms)
 - Depth 8.8-inches (22.4cms)
 - Height 13.9-inches (35.3cms)

CAPACITY

- Operator 500lbs (227kgs)
- Basket 25lbs (11.4kgs) / 1933 cubic inches
- Total Capacity 525lbs (238kgs)
- Load-tested 1200+lbs (544.3+kgs)

SAFETY

- Battery overcharge protection
- Electromagnetic braking w/manual release
- Dynamic regenerative braking
- Automatic parking brake
- Rear safety bumper
- Front and rear stability devices
- Drive-off charging protection
- ASIC micro-processor controller
- Flip-up/adjustable/removable arms (optional)

WARRANTY**

- Platform, seat post – Lifetime
- Transaxle, motor, controller, battery charger – 2 years
- Seat – 1 year
- Batteries, labor – 6 months

PROVIDES

- Ergonomically designed handle
- Protected throttle operation
- ON/OFF Switch
- All steel one-piece extended platform/frame, powder coated
- Solid, flat-free, non-marking tires
- Front tire 6-inches (16.2cms)
- Rear tires 8-inches (20.3cms)
- Custom colors available (SmartMover only)

**Optimum range. Several factors affect range such as payload, battery type and age, charging habits, terrain traveled, ambient temperature, etc. **Refer to actual warranty for complete details. Specifications subject to change without notice.*

AMIGO MOBILITY INTERNATIONAL, INC.
6693 Dixie Highway • Bridgeport, MI 48722-9725
800 248-9131 • www.myamigo.com