

Increase Customer Satisfaction with Electric Convenience Vehicles

Provide safe and convenient transportation for visitors and employees

Increase: *Customer satisfaction, staff productivity*

Reduce: *Slip and falls, labor costs*

Attract: *New visitor usage, community awareness*



"We've heard nothing but positive comments about our new Amigos from those who use them. The Amigos provide our visitors and guests with unprecedented freedom of movement. Now, instead of being dependent on other family members or staff to get around, our guests direct their own travel throughout our facility. The Amigos are reliable, maneuverable and easy to use. Our only complaint is that we didn't order them sooner."

*University of Michigan Hospital
Security and Entrance Services*



Providing Electric Convenience Vehicles is a Smart Investment

Amigo Mobility Offers 40 Years of "Improving Lives Through Mobility"

Amigo Features

- All-Through Steel Platform Resists Cracks/Breaks
- Long-term Commercial Durability Medical Equipment with More Plastic/Less Steel
- 49" Platform Length Fits Broad Range of People
- 4" Longer Platform Option (SmartMover XL)
- Virtually Indestructible Seat
- Permanently Attached Commercial Basket
- 500 lbs. Rider Weight Capacity

Safety Features

- Programmable Speed Control
- Protected Throttle Operation
- Optional Seat Switch
- Battery Overcharge Protection

Optional Features

Four Basket Selections: Trimline, Mid-Size, Large or Locked Basket with Cover (*good for transporting medication*)

Dimensions	Trimline	Mid-Size	Large
basket width at top	16	21.25	25.5
basket width at bottom	15	18.75	21
basket length at top	9	15.25	20
basket length at bottom	8	12.75	16.5
basket height at front	14.5	9.25	15.25
basket height at back	14.5	9.25	11
Cubic Inches	1,600	2,611	6,200

Locked Basket with Cover

Height 20
Width 28"
Depth 14"



- Arms
- Front Manual Horn/Reverse Alert Horn
- Cane Holder
- O2 Tank Holder
- Swivel Seat w/Locking Mechanism
- Customizable Décor and Color Scheme

AMIGO MOBILITY INTERNATIONAL, INC.
1-800-248-9130 • www.myamigo.com

© 2007 Amigo Mobility International, Inc. Amigo® is a trademark of Amigo Mobility International, Inc.

Like No Other Product...

- Made in the USA
- Completely Recyclable
- Low Cost of Ownership
- Responsive Technical Service with Live Support
- Limited Lifetime Warranty



Your Local Representative

M800-012108