



# TRAVELMATE RD

## SPECIFICATIONS MODEL 820000

### TRAVEL SURFACES

- Grass, gravel, sidewalks, linoleum, carpeting, indoor/outdoor surfaces

### SPEED AND POWER

- 0-4.2 mph (0 - 6.8 kph)
- 16 miles (25.7 km) battery capacity\*
- 24V electrical system for two 12AH, airline-safe, sealed batteries
- 2A off board battery charger w/AC Cord

### UNIPAC® & VARI-SPEED DRIVE SYSTEM

- Sealed, Lifetime-Lubricated, Rear Drive
- No Belts / No Chains

### MANEUVERABILITY

- 35 inch (89 cm) turning radius
- 2 inch (5 cm) underbody ground clearance

### DIMENSIONS

- 45 inches (114.3 cm) long
- 22.5 inches (57.2 cm) rear wheel width
- 36 inches (91.4 cm) high
- 94 pounds (42.6 kg) (with batteries)

### COLOR

- Moss Green

### SAFETY

- MPC Electronics Controller
- Electromagnetic Braking with Manual Release
- Seat Lock
- Flip-Up, Adjustable, Removable Arms (Optional)
- Rear Safety Bumper with Anti-Tip Wheels
- 26 inch (66 cm) infinitely adjustable EZ-Tilt Handle with HI-LO Controls

### PROVIDES

- Steel, powder coat, one piece platform
- Dynamic, regenerative & electromagnetic brakes
- 8" x 2" (20.3 x 5.1 cm) flat-free front tire
- 9" x 3" (22.9 x 7.6 cm) flat-free rear tires
- 360-degree seat rotation
- 250 pounds (113.4 kg) weight capacity
- Programmable Controller
- Amigo Limited Lifetime Warranty

*\*Range based on a 180-pound (81.6 kg) driver on a flat surface. Range depends on battery type, driver weight, terrain traveled, ambient temperature, and user's charging routine.*

*Products with options as shown are subject to change without notice.*

**AMIGO MOBILITY INTERNATIONAL, INC.**  
6693 Dixie Highway • Bridgeport, MI 48722-9725  
800 248-9131 • [www.myamigo.com](http://www.myamigo.com)

M1364-09